## Onse Step Morphine MonlabTest®

IVD

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

### INTENDED USE

The One-Step Morphine MonlabTest® is a lateral flow, one-step immunoassay for qualitative detection of morphine and their metabolites in human urine at a cut-off of 300ng/mL. This product is used to obtain a visual, qualitative result and is intended for professional use. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary results are indicated.

### SUMMARY

The opiates such as heroin, morphine, and codeine are derived from the resin of opium poppy. Heroin is quickly metabolized to morphine. Thus, morphine and morphine glucuronide both might be found in the urine of a person who has taken only heroin. The body also changes codeine to morphine. Thus, the presence of morphine (or the metabolite, morphine glucuronide) in the urine indicates heroin, morphine and/or codeine use.

Urine based screening tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The One-Step Morphine MonlabTest® is based on the principle of highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids. This test is a rapid, visual, competitive immunoassay that can be used for the qualitative detection of morphine in human urine at a 300ng/mL cut-off concentration.

### PRINCIPLE

The One-Step Morphine MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug-protein conjugate) competes with the drug which may be present in urine for limited antibody binding sites. The test device contains a membrane strip which is precoated with drug-protein conjugate on the test band. The colored anti-morphine antibody-colloidal gold conjugate pad is placed at the end of the membrane. In the absence of drug in the urine, the solution of colored antibody-colloidal gold conjugate and urine moves upward chromatographically by capillary action across the membrane. This solution then migrates to the immobilized drug-protein conjugate zone to form a visible line as the antibody complexes with the drug-protein conjugate. Therefore, the formation of a visible precipitant in the test band region occurs when the test urine is **negative** for the drug.

When the drug is present in the urine, the drug/metabolite antigen competes with the drug-protein conjugate on the test band region (T) for the limited antibody sites on the anti-morphine antibody-colloidal gold conjugate.

When a sufficient concentration of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug-protein conjugate zone on the test band region. Therefore, absence of the color band on the test region indicates a **positive** result.

A control band with a different antigen/antibody reaction is also added to the immunochromatographic membrane strip at the control region (C) to indicate that the test is performed properly. This control line should always appear regardless of the presence of drug or metabolite. This means that **negative** urine will produce **two** colored bands, and **positive** urine will produce only **one** band. The presence of this colored band in the control region also serves as verification that 1) sufficient specimen has been added, and 2) that proper flow was obtained.

### REAGENTS AND MATERIALS SUPPLIED

- 40 individually wrapped test devices which include one disposable pipette each.
- One instruction sheet.

### MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer

### STORAGE AND STABILITY

The test kit should be stored at temperature 2-30°C. Each device should remain in its sealed pouch for the duration of the shelf life of two years.

### PRECAUTIONS

- FOR IN-VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Do not use it after the expiration date.


### SPECIMEN COLLECTION AND HANDLING

The One-Step Morphine MonlabTest® is formulated for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing may be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing. Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.

### TEST PROCEDURE


- Review "Specimen collection" instructions. Test device, samples of patient, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control number.
- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15ml) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- **Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

## INTERPRETATION OF RESULTS


### Negative

Two pink-rose lines (bands) are visible in the control (C) and test (T) areas of the test window. **The intensity of the test line may be less than that of the control line; this still means negative result.**

### Positive

The control line appears in the test window, but the test line is not visible.

### Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

**Note:** A very faint line in the test region indicates that the THC in the samples is near to the cut-off level for the test. These should be re-tested or confirmed with a more specific method before a positive determination is made.

## LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical and/or procedural errors as well as other substances or factors not listed may interfere with the test and cause false results. See **SPECIFICITY** for lists of substances that will produce positive result, or that do not interfere with test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

## QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance.

Before using a new kit with patient specimens, positive and negative controls should be tested. Quality control specimens are available from commercial sources. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

## PERFORMANCE CHARACTERISTICS

### Accuracy

The accuracy of the One-Step Morphine MonlabTest® was evaluated in comparison to other commercially available methods. One hundred fifty (150) urine samples, collected from presumed non-user volunteers, have been tested by both tests; all were found negative by both methods with 100% agreement.

In a separate study, fifty (50) urine samples, obtained from a clinical laboratory and determined morphine positive by GC/MS, were tested with One-Step Morphine MonlabTest® in comparison to other commercially available methods.

Concentration (ng/mL)	Number Tested	commercially available kit (300ng/mL)		MonlabTest® (300ng/mL)	
		Pos.	Neg.	Pos.	Neg.
0–250	25	0	25	0	25
320–1500	25	25	0	25	0

### Reproducibility

The reproducibility of the One-Step Morphine MonlabTest® was evaluated at four different sites using controls. Of the fifty (50) samples with morphine a concentration of 150ng/mL, all were determined negative. Of the fifty samples with a morphine concentration of 600ng/mL, all were determining positive.

### Precision

The precision of the One-Step Morphine MonlabTest® was determined by conducting the test with spiked controls. The control at 150ng/mL gave negative results and the control at 600ng/mL gave positive results.

### Specificity

The specificity for the One-Step Morphine MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

**The following structurally related compounds produced positive results when tested at levels equal to or greater than the concentrations listed below.**

Compound	Concentration (ng/mL)
Morphine	300
Codeine	300
Ethyl Morphine	300
Hydrocodone	375
Hydromorphone	400
Morphine-3-β-d-glucuronide	490

**The following compounds were found not to cross-react when tested at concentrations up to 100µg/mL.**

Acetaminophen	Imipramine
Cyclobenzaprine	Phencyclidine
Acetone	(+/-)-Isoproterenol
r-Cyclodextrin	Phenelzine
Albumin	Lidocaine
C-yproheptadine	Pheniramine
Amitriptyline	Maprotiline
(-)-Desoxyephedrine	Phenobarbital
D-Amphetamine	Mephentermine
Dextromethorphan	Penothiazine
L-Amphetamine	MDA
Diazepam	Meperidine
Amoxaoine	Mephentermine
4-Dimethylaminoantipyrine	β-Phenylethylamine
Ampicillin	Mes
Diphenhydramine	Phentermine
Aspartame	Methadol
(Asp-Pheethyl-Ester)	Primidone
5,5-Diphenylhydantoin	Methadone
Aspirin	Procaine
(Acetylsalicylic Acid)	Methamphetamine
Dopamine	Promazine
(3-Hydroxytyramine)	Methapyrilene
Atropine	Promethazine
Doxylamine	Methaqualone
Baclofen	d-Propoxyphene

Ecgonine	(1R,2S)-(-)-N-Methyl-ephedrine
Ecgonine Methyl Ester	Protriptyline
Benzocaine	Methylphenidate
(Ethyl-Aminobenzoate)	Pseudoephedrine HCL
(-)-Ephedrine	(+/-)-3,4-Methylenedioxy-
Benzoylcegonine	2-propylantanoic acid
(+/-)-Ephedrine	methamphetamine(+/-)-3,4-MD
Bilirubin	MA Quinidine
(-)-ψ-Ephedrine	Naloxone
(+)-Brompheniramine	Quinine
(Dexbrompheniramine)	Naltrexone
(+)-ψ-Ephedrine	Ranitidine
Caffeine	β-Naphthaleneacetic Acid
(+/-)-Ephedrine	Riboflavin
Carbamate	(+)-Naproxen
Erythromycin	Salicylic acid
Carbamazepine	(-)-Nicotine
Ethanol	Secobarbital
Carisoprodol	Nicotinic Acid
Chloroquine	Sodium Chloride
(+)-Chlorpheniramine	(+/-)-Norephedrine
Fentanyl	((+/-)-Phenylpropanol)
(+/-)-Chlorpheniramine	Sulindac
Fenofibrate	Nortriptyline
Chlorprothixene	Theophylline
Gemfibrozil	Noscapine
Chlorthalidone	cis-Thiothixene
Glucose	Oxalic Acid
Cocaine	Trifluoperazine
Guaiaicol glyceryl ether	Oxycodone
Clofibrate	Trimethobenzamide
DL-Homatropine	Penicillin-G(Benzylpenicillin)
Creatine	Tyramine
Hemoglobin	Pentobarbital
Creatinine	Vitamin C
Ibuprofen	Perphenazine

## REFERENCES



- Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn. M. J. And Barce1oux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews . Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B,1987

## PACKAGING

CHEM-7001C-40

40 Morphine Tests

## GRAPHICAL SYMBOLS USED

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

**Manufacturer:** Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

 **European Representative:**

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

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## REFERENCIAS

- Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn. M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

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**Fabricante:** Shanghai Chemtron Biotech Co., Ltd.