NOROVIRUS MonlabTest® MO-076004 20 TESTS

One-step immunochromatographic test for the differential detection of
Norovirus Genogroups I and II in human faeces

For *in vitro* use only. Store at 2° to 30°C.

PRINCIPLE

The Norovirus MonlabTest is a chromatographic immunoassay for the qualitative detection, in separate bands, of Genogroup I (GI) and Genogroup II (GII) of Norovirus in human faeces. A positive signal in either test bands provides a good indication that we may be in presence of a Norovirus infection which should draw the attention of the clinician. The test is based on the immunological capture of coloured microbeads during its moving along a membrane on which specific monoclonal antibodies against GI and GII have been immobilized in two separate bands.

Norovirus MonlabTest uses a combination of:

- 1) Red latex particles conjugated to specific antibodies against GII which cooperates with antibodies specific for GII located on the membrane under the control band.
- 2) Other red latex particles conjugated to specific antibodies against GI which cooperates with antibodies specific for GI located on the membrane above the control band.
- 3) Blue latex particles conjugated to an antigen recognized by an antibody specific for this antigen bound to the membrane forming the so called test control band.

In this test the sample is treated first with the sample diluent buffer (included in the kit) to achieve the extraction of the virus from the stool matrix. After the extraction, you just need to add a certain volume of supernatant onto the test strip and wait for 15 minutes.

When the extracted sample flows through the test membrane, the coloured particles migrate. In the case of a positive sample, specific antibodies on the membrane will capture the coloured particles covered by the antigen.

Depending on the virus content of the sample, different coloured lines will be visible. These lines are used for interpretation of results, following 15 minutes incubation at room temperature (see Figure 1).

INTENDED USE

DESCRIPTION OF NOROVIRUS

Noroviruses are a kind of single-stranded, positive-sense RNA viruses belonging to the family *Caliciviridae*^{1,2,3} (Sapoviruses also belong to this family, among others). They are highly contagious viruses whose main routes of transmission are: person-to-person contact and through contaminated food/water (in the U.S. is estimated that they account for 50% of gastrointestinal outbreaks by food poisoning). They often cause large outbreaks in closed communities in a variety of environments such as hospitals, nursing homes, schools, kindergarten, restaurants, cruises, ... where, once the virus has been introduced, the infection spreads very quickly. Several studies show that Noroviruses are responsible for almost 50% of gastroenteritis outbreaks worldwide (considering all outbreaks of any aetiology)³.

Human Noroviruses have been very difficult to study because of their high diversity and the lack of an efficient cell culture for *in vitro* replication as well as a suitable animal model. The cloning of genomes of Norovirus and expression of the viral capsid proteins in baculovirus and other expression systems, have allowed to obtain "virus-like particles" (VLP) which, once assembled, are able to recreate the real structure of the virus.

Norovirus capsid is composed of a single major structural protein, the capsid protein (VP1) which can be divided into two main domains: the shell (S) and the protruding domains (P)⁴. The expression of VP1 in a eukaryotic system allows to obtain empty structural particles similar to the native present in the virus (VLP) that have been used as a substitute for native Norovirus for many years in any research process⁵.

CLASSIFICATION OF NOROVIRUSES

Noroviruses are grouped into five Genogroups (GI to GV), of which GI and GII are involved in most acute cases of viral gastroenteritis in humans.

Within each Genogroup viruses are classified into Genotypes (GI.1, GII.1...). Over 25 different Genotypes have been described within Genogroups I and II. Of these, GII.4 is the most common Genotype representing around 60-80% of cases worldwide^{6,7}. Following are GII.2, GII.3 and GII.7 within GII and GI.1, GI.3 or GI.4 within GI⁴ although, in general, this Genogroup is rather less common than GII, representing less than 5%⁸.

Recent studies have shown that in recent years new variant strains of GII.4 have been identified, which arise as a result of genetic changes that, in some cases, affect only one amino acid in the sequence of the capsid protein but this is enough to make them different from previous GII.4 strains. This means that any type of immunity that an individual would have generated against GII.4 could be useless against these new strains. In short, a single Genotype is presented as a whole family, tremendously complex and variable over time⁷.

CHARACTERISTICS OF NOROVIRUS INFECTION

Although Noroviruses can be detected throughout the year, it has been clearly observed that there is a seasonal prevalence with peaks during autumn and winter³.

Norovirus infection has an incubation period of 24-48 h and is characterized by nausea, vomiting, abdominal pain, fever, ... Acute symptoms usually appear for 1-3 days, although safety hygiene actions must be maintained for a minimum of two weeks (the elimination of the virus in faeces may continue for days, even weeks); the patient may remain asymptomatic from third or fourth day after onset of symptoms⁹.

After overcoming the infection, immunity to Norovirus is usually incomplete and temporary (about 6 months) as well as specific for a particular Genotype¹⁰. Given the high genetic variability of Noroviruses, individuals are likely to be infected several times throughout their lifetime. This may explain the high rates of infection that occur in all ages at an outbreak level. Recent studies suggest that susceptibility to Norovirus infection may be genetically determined¹¹.

MATERIALS PROVIDED

- 20 cassettes
- 20 vials with diluent (1,5mL)
- 20 non-graduated disposable plastic pipettes (yellow)
- Instructions for use

MATERIALS NOT INCLUDED IN THE KIT

- Vortex
- Stopwatch

PRECAUTIONS

1. Patient specimens (faeces) should be handled with caution as they may contain infectious agents. All the necessary protections should be used throughout handling.
2. The sample dilution buffer contains sodium azide as anti-microbial agent. Avoid direct contact with skin and mucous membranes. Dispose of appropriately. Do not use the buffer if signs of contamination are present.
3. Do not eat, drink, smoke or prepare or store food in areas where reagents or samples are handled.
4. Once completed the procedure, clean work surfaces with soap and water and finish disinfecting with a suitable solution. Finally, discard the gloves and then wash hands first with soap and water rubbing them well.
5. Do not exchange components from one kit lot to another.
6. Before use, allow all kit components and stool samples to reach room temperature, as cold reagents and/or samples can reduce test functionality. 20 to 30 minutes are recommended to reach room temperature.
7. Use the reagents only *in vitro*.
8. Do not use kit components beyond their expiration date.
9. In case of package damage, the product may be used if none of its components has been damaged.
10. It's very important to add the correct volume of extracted sample to the reaction device. If the amount is lower than required, the chromatography may not be completed as the sample does not reach the reaction area, if it is bigger, brown lines may appear instead of red or blue lines.
11. Used product should be disposed of according to applicable laws.
12. Do not use the test if any coloured lines appear in the result area before using the test.
13. It is critical to collect the correct sample quantity: approximately **110 mg** of a **solid sample** (a **small ball of 5 mm in diameter**). If the **sample is semi-liquid** (unable to take it with a pipette), take a sample amount so that it completely **covers the grooves of the stick** attached to the vial cap and **110 µl** of a liquid sample (**4 drops** if the disposable pipettes included in the kit are used). These quantities are extracted in the sample diluent supplied with the vials provided in the kit. An excess of sample in relation to the amount of buffer added prevents the chromatography from running correctly; this is especially critical in the case of solid samples, since it is harder to obtain an appropriate quantity.
14. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and lotification.

STORAGE AND STABILITY

Product may be stored at any temperature between 2 and 30°C. The expiry date is printed on the tubes or aluminium packages.

SAMPLES

- This test is designed for the analysis of human faecal samples.
- It is recommended to collect the stool sample as soon as symptoms appear (diarrhoea and vomiting, especially) as the elimination of virus in faeces is greatest during the first three days from the onset of symptoms.
- Do not use samples collected in transport media or added with preservatives (such as formalin, SAF, PVA or similar) or enrichment media as their presence could interfere with the proper performance of the test.
- The use of **fresh untreated samples** is highly recommended. If they have to be stored for some time, they can be stored in the refrigerator (+2-8°C) for 1 or 2 days. For a longer storage, they should be frozen at -20°C, considering that some samples lost immunoreactivity after being frozen.
- Pay special attention when analyzing haemorrhagic samples as they often give false positive results when the blood content is high. An indicator of this destabilization of the test is usually the alteration of the colour in the control band (tends to show a purple or dark blue colour).
- In case of freezing, thaw the samples completely at room temperature prior to analysis. Avoid freezing-thawing cycles for the faecal samples as the immunological recognition of virus can be altered.

PREPARATION OF FAECAL SAMPLES

General Note: all the necessary protections should be used throughout the test procedure due the handling of infectious samples. Once the work is concluded, do not forget to comply with the hygiene procedures detailed in point 4 of the "Precautions" section.

1. Homogenize previously the sample to get an aliquot as much representative as possible.
2. Unscrew the cap from the vial with caution in order not to spill the sample diluent buffer. With **solid samples**, take with the stick attached to the vial cap an approximate amount of **110 mg** of faeces (a small portion of **5 mm in diameter**). If the sample is **semi-liquid** (unable to take it with a pipette), take a sample amount so that it completely **covers the grooves of the stick** attached to the vial cap. With **liquid samples**, take a volume of **110 µl** (**4 drops** if the disposable pipettes included in the kit are used).
3. Carefully add the sample into the vial containing the dilution buffer. Screw the cap well and shake vigorously to ensure a homogeneous mixture.

PROCEDURE

Once the samples have been prepared as described in previous section, the procedure is as follows:

1. Remove the reaction device from the aluminium pouch. Discard the desiccant bag as it only serves to preserve the test from moisture.
2. Break the top end of the vial stopper.
3. Put the vial upside down and add **3 drops** to the sample area of the reaction device (round window marked with an arrow).
4. Wait for **15 minutes**, read and interpret the results.

RESULTS

The five results shown in Figure 1 are an example of the results that can be obtained with the Norovirus MonlabTest.

Three different coloured bands can be distinguished:

Blue band: is the control band indicating a proper test performance.

Upper red band: indicates the presence of Norovirus Genogroup I in the sample.

Lower red band: indicates the presence of Genogroup II in the sample.

The control blue band should always appear. The additional presence of any red band indicates the presence of Norovirus in the sample.

Result 1: NEGATIVE result: a single transverse **BLUE** line appears in the middle area of reaction device along with the letter "C" marked on the plastic cassette. This band is the control band and should always appear indicating that chromatography takes place normally.

Result 2-4: POSITIVE results

Result 2. Detection of GII: a **BLUE** band (control) and a **RED** band 3 mm below the control band aligned with the label "T1" marked on the casing appear. The band intensity depends on virus concentration in the sample.

Result 3. Detection of GI: a **BLUE** band (control) and a **RED** band 3 mm above the control band aligned with the label "T2" marked on the casing appear. The band intensity depends on virus concentration in the sample.

Result 4. Detection of GI and GII: a **BLUE** band (control) along with two **RED** bands (one above -GI- and one below -GII- the control band) appear.

Result 5: INVALID results: the blue control band does not appear or the blue colour of this band is clearly altered (very dark blue or a purple colour), also nonspecific colours in the positive bands (different from red) are invalid results. In general, any colour combination different from those indicated in results 1-4 indicates an anomalous test performance.

Possible reasons are:

- Some reagents have got damaged or the test has expired.
- The sample was not prepared according to the instructions of use.
- High blood content in the sample.

In the event of an invalid result it is recommended that another test is run, strictly following the protocol described in this manual. For bloody samples, it is advisable the use of an alternative technique because the problem of destabilization does not usually depend on the cassette but of the sample matrix.


Figure 1: Examples of the potential results.

Any line that appears beyond 15 minutes because of the nature of the sample has no diagnostic value.

NOTE: the final and definitive diagnosis is established by the clinician. This test only detects Norovirus in a sample but is not an argument to conclude that the person has an infection with Norovirus.

WARNING: Including our controls with an established result is recommended to ensure that the data obtained are correct.

PROCEDURE LIMITATIONS

1. The Norovirus MonlabTest is intended for the differential identification of Norovirus Genogroups I and II by detecting its presence in human faeces provided that the viral load is equal to or higher than the detection limit of the assay.
2. This is a qualitative, not quantitative test, although the intensity of positive bands relates to the amount of detectable virus in the stool sample.
3. Over 200 stool samples were assessed to ensure proper test performance. The correlation of results with other techniques (RT-PCR, ELISA and rapid tests) was good. However, this study does not exclude potential interference in test performance when analysing other faecal samples.

4. Norovirus MonlabTest has not been validated with all Genotypes able to infect humans and, therefore, the test may fail to detect some Norovirus due to the high antigenic diversity of the circulating strains.
5. Insufficient sample amount may lead to extremely weak positive results. In this case, repeat the test with a larger amount of sample while maintaining the recommended ratio to the volume of sample diluent. On the other hand, a sample excess can cause the test to run very slowly and even prevent the test to perform properly (the control line is not visible). In this event, the test must be repeated with a lower sample amount (see section "Preparation of Faecal Samples").
6. A negative result does not exclude the possibility of a Norovirus infection. The failure in the detection of Norovirus may be the result of factors such as: taking the sample at an incorrect time of disease (when a low amount of virus is eliminated in the faeces), an incorrect sample storage, an inadequate handling of the sample, the presence of a Norovirus Genotype not detected by the strip (as the test has not been validated with all Genotypes able to infect humans and may fail to detect any of them).
7. A positive result does not exclude the presence of other pathogens; furthermore, a co-infection of Norovirus with other microorganisms may occur. In this sense, note that the most common co-infections of Norovirus often occur with parasites, especially with *Giardia lamblia* followed by *Cryptosporidium parvum*¹². In any case, co-infections can only be clarified by a differential diagnosis.
8. It is observed that faecal samples with high blood content may interfere negatively with the test. In this case, specificity problems may appear with Norovirus negative samples. This destabilization of the test is often accompanied by an alteration in the colour of the control band; a purple or dark blue colour appears instead of the expected light blue (see images in "Results Reading" section).
9. Test results should be interpreted in conjunction with information available from epidemiological studies, patient clinical evaluation and other diagnostic procedures.

SENSIBILITY AND SPECIFICITY

ANALYTICAL SENSITIVITY

The internal standard to validate the manufactured lots is a mixture of "virus like particles" GI.1+GII.4 as they are the most representative and common genotypes within each genogroup.

A mean sensitivity of 6.25 ng/mL is obtained for NoV GI.1 and a mean sensitivity of 0.75 ng/mL for NoV GII.4 though up to 1.5 ng/mL and 0.1 ng/mL are usually detected for NoV GI and NoV GII respectively.

It is important to note that the detection limit of the test was also analyzed using real stool samples. Values were consistent with those obtained from internal standards. These results demonstrate the robustness of the test.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The sensitivity and specificity of the Norovirus MonlabTest was assessed by measuring a total of **96 negative samples**, **23 positive GI samples** and **100 positive GII samples**, characterized by RT-PCR (reference technique). All samples analyzed were frozen.

The following results were obtained:

Total Sensitivity: 85%

GI Sensitivity: 87%

GII Sensitivity: 85%

Specificity: 96%

It can be seen that the specificity of the Norovirus test is very high, higher than 95%. With respect to sensitivity, it is very important to analyse the values with caution taking into account that the reference technique (RT-PCR) is ultra-sensitive and sometimes it might detect as positive samples belonging to asymptomatic patients or patients who have overcome the disease. Anyway, sensitivity values around 85% for both GI and GII, are high ones for a rapid test.

REPEATABILITY

Ten replicates of each of the three concentrations established as PC ("positive control"), LPC ("low positive control") and NC ("negative control") with our internal standard were measured on the same day by the same person. A 100% repeatability was obtained with these three critical concentrations indicating a high intra-assay precision of the test.

REPRODUCIBILITY

INTER-DAY PRECISION

Using a single lot of Norovirus test a sensitivity curve is measured through four days spaced in time. The results were very reproducible (the same sensitivity for both GI and GII through the four days of measurement).

INTER-OPERATOR PRECISION

Five operators measured in duplicate a sensitivity curve. Differences were observed but in no case exceeded 1 two-fold dilution.

INTER-LOT PRECISION

Using three different lots of Norovirus test a sensitivity curve was measured in duplicate. The analysis was performed by the same person on the same day. Only differences below to 1 two-fold dilution were observed, acceptable and tolerable for the assay.

The differences found in the different "Reproducibility" sections are acceptable for a qualitative immunochromatographic technique with its inherent variability.

HOOK EFFECT

The maximum amount of Norovirus that a person can eliminate during the acute phase of disease is about 10^{12} particles/g of stool¹³, which is equivalent to 7.5×10^{10} particles/mL taking into account that Norovirus MonlabTest uses a sample extraction of 75 mg of stool in 1 mL of buffer. Considering an approximate equivalence between VLPs and mass of: 1 µg equal to 5.87×10^{10} VLPs/mL, 7.5×10^{10} particles/mL are about 1280 ng/mL of VLP. For excluding Hook Effect, a maximum concentration of 40000 ng/mL of both GI.1 and GI.4 was measured with the Norovirus test and, in both cases, positive results were obtained. These concentrations are about 6400-fold the mean sensitivity reached for GI.1 and about 53300-fold the mean sensitivity reached for GI.4 (see section "Analytical Sensitivity"). It is concluded that Hook Effect is not a problem for this test.

INTERFERING SUBSTANCES

The substances indicated in the below table at the concentration specified did not interfere with the results of the test when added to stool samples (positive and negative ones):

Racecadotril	5% (p/v)	Ibuprofen	20% (p/v)
Cimetidine	10% (p/v)	Acetylsalicylic acid	30% (p/v)
Loperamide	5% (p/v)	Edulcorant	5% (p/v)
Metronidazole	10% (p/v)	Palmitic acid	40% (p/v)
Omeprazole	3% (p/v)	Barium Sulphate	5% (p/v)
Ampicillin	15% (p/v)	Mucin	5% (p/v)

CROSS REACTIVITY

The Norovirus MonlabTest was evaluated against different micro-organisms likely to be present in the intestinal tract at any time in a sufficiently high concentration. No cross-reactivity was detected with any of the following microorganisms:

BACTERIA:

Aeromonas baumannii, *Aeromonas hydrophila*, *Aeromonas caviae*, *Bacillus spp.*, *Burkholderia cepacia*, *Campylobacter coli*, *Campylobacter jejuni*, *Citrobacter freundii*, *Clostridium perfringens*, *Clostridium difficile*, *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Enterococcus faecium*, *Escherichia coli*, *Escherichia coli O111*, *Escherichia coli O127*, *Escherichia coli O26*, *Escherichia coli O55*, *Escherichia coli O157:H7*, *Hafnia alvei*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Lactobacillus casei BL23*, *Lactococcus lactis MC1363*, *Listeria monocitogenes*, *Morganella morganii*, *Plesiomonas shigelloides*, *Proteus mirabilis*, *Proteus vulgaris*, *Proteus penneri*, *Providencia stuartii*, *Pseudomonas aeruginosa*, *Pseudomonas stutzeri*, *Salmonella choleraesuis*, *Salmonella enteric serogroup B*, *Salmonella enteric serogroup D*, *Salmonella typhi*, *Serratia marcescens*, *Shigella flexnerii*, *Shigella sonnei*, *Staphylococcus aureus*, *Streptococcus agalactiae*, *Streptococcus dysgalactiae*, *Streptococcus viridans*, *Vibrio cholerae*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*.

VIRUSES:

Astrovirus, *Adenovirus*, *Enterovirus*, *Rotavirus strain Wa*, *Rotavirus*, *Sapovirus*, *virus Aichi*.

FUNGI/PARASITES/OTHER:


Blastocystis hominis, *Candida albicans*, *Candida parapsilosis*, *Candida tropicalis*, *Cryptosporidium parvum*, *Entamoeba histolytica*, *Entamoeba coli*, *Giardia lamblia*.

REFERENCES

- Atmar RL and Estes MK. *Diagnosis of non-cultivable gastroenteritis viruses, the human caliciviruses*. Clinical Microbiology Reviews, 2011. Vol: 14 (pp: 15-37).
- Ribes Fernández JM and Buesa Gómez J. *Infecciones por Norovirus*. Enferm Infecc Microbiol Clin, 2010. Vol: 28 (pp: 51-55).
- Patel MM, Hall AJ, Vinjé J and Parashar UD. *Noroviruses: a comprehensive review*. Journal of Clinical Virology, 2009. Vol: 44 (pp: 1-8).
- Xerry J, Gallimore CI, Iturriza-Gómara M and Gray JJ. *Genetic characterization of Genogroup I Norovirus in outbreaks of gastroenteritis*. Journal of Clinical Microbiology, 2010. Vol: 48 (pp: 2560-2562).
- Ausar SF, Foubert TR, Hudson MH, Vedvick TS and Middaugh CR. *Conformational stability and disassembly of Norwalk virus-like particles*. The Journal of Biological Chemistry, 2006. Vol: 281 (pp: 19478-19488).
- Buesa J, Montava R, Abu-Mallouh R, Fos M, Ribes JM, Bartolomé R, Vanaclocha H, Torner N and Domínguez A. *Sequential evolution of Genotype GI.4 Norovirus variants causing gastroenteritis outbreaks from 2001 to 2006 in Eastern Spain*. Journal of Medical Virology, 2008. Vol: 50 (pp: 1288-1295).
- Siebenga J, Kroneman A, Vennema H, Duizer E and Koopmans M. *Food-borne viruses in Europe network report: the norovirus GI.4 2006b (for US named Minerva-like, for Japan Kobe034-like, for UK V6) variant now dominant in early seasonal surveillance*. Eurosurveillance, Jan-Mar 2008. Vol: 13 (Issues 1-3).
- Okame M, Akihara S, Hansman G, Hainian Y, Tran HT, GiaPhan T, Yagyu F, Okitsu S and Ushijima H. *Existence of multiple Genotypes associated with acute gastroenteritis during 6-year survey of Norovirus infection in Japan*. Journal of Medical Virology, 2006. Vol: 78 (pp: 1318-1324).
- Furuya D, Kuribayashi K, Hosono Y, Tsuji N, Furuya M, Miyazaki K and Watanabe N. *Age, viral copy number, and immunosuppressive therapy affect the duration of norovirus RNA excretion in inpatients diagnosed with norovirus infection*. Jpn. J. Infect. Dis, 2011. Vol: 64 (pp: 104-108).
- Reeck A, Kavanagh O, Estes MK, Opekun AR, Gilger MA, Graham DY and Atmar RL. *Serological correlate of protection against Norovirus-induced gastroenteritis*. Journal Infection Disease, 2010. Vol: 202 (pp: 1212-1218).
- Waqas Nasir. *A study of Norovirus-HBGA interactions*. Chalmers University of Technology. University of Gothenburg, Sweden. December 2009.

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SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufactured by		For in vitro diagnostic use
	Do not re-use		Please read pack insert
	Contains sufficient for <n> tests		Dry storage
	Catalogue number		Store at
	Lot number		Expiry date
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices		Dilution buffer
	Reagent		

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


- Atmar RL and Estes MK. *Diagnosis of non-cultivable gastroenteritis viruses, the human caliciviruses*. Clinical Microbiology Reviews, 2011. Vol: 14 (p. 15-37).
- Ribes Fernández JM and Buesa Gómez J. *Infecciones por Norovirus*. Enferm Infecc Microbiol Clin, 2010. Vol: 28 (pp: 51-55).
- Patel MM, Hall AJ, Vinjé J and Parashar UD. *Noroviruses: a comprehensive review*. Journal of Clinical Virology, 2009. Vol: 44 (pp: 1-8).
- Xerry J, Gallimore CI, Iturriza-Gómara M and Gray JJ. *Genetic characterization of Genogroup I Norovirus in outbreaks of gastroenteritis*. Journal of Clinical Microbiology, 2010. Vol: 48 (pp: 2560-2562).
- Ausar SF, Foubert TR, Hudson MH, Vedvick TS and Middaugh CR. *Conformational stability and disassembly of Norwalk virus-like particles*. The Journal of Biological Chemistry, 2006. Vol: 281 (pp: 19478-19488).
- Buesa J, Montava R, Abu-Mallouh R, Fos M, Ribes JM, Bartolomé R, Vanaclocha H, Torner N and Domínguez A. *Sequential evolution of Genotype GII.4 Norovirus variants causing gastroenteritis outbreaks from 2001 to 2006 in Eastern Spain*. Journal of Medical Virology, 2008. Vol: 50 (pp: 1288-1295).
- Siebenga J, Kroneman A, Vennema H, Duizer E and Koopmans M. *Food-borne viruses in Europe network report: the norovirus GII.4 2006b (for US named Minerva-like, for Japan Kobe034-like, for UK V6) variant now dominant in early seasonal surveillance*. Eurosurveillance, Jan-Mar 2008. Vol: 13 (Issues 1-3).
- Okame M, Akihara S, Hansman G, Hainian Y, Tran HT, GiaPhan T, Yagyu F, Okitsu S and Ushijima H. *Existence of multiple Genotypes associated with acute gastroenteritis during 6-year survey of Norovirus infection in Japan*. Journal of Medical Virology, 2006. Vol: 78 (pp: 1318-1324).
- Furuya D, Kuriyayashi K, Hosono Y, Tsuji N, Furuya M, Miyazaki K and Watanabe N. *Age, viral copy number, and immunosuppressive therapy affect the duration of norovirus RNA excretion in inpatients diagnosed with norovirus infection*. Jpn. J. Infect. Dis, 2011. Vol: 64 (pp: 104-108).
- Reeck A, Kavanagh O, Estes MK, Opekun AR, Gilger MA, Graham DY and Atmar RL. *Serological correlate of protection against Norovirus-induced gastroenteritis*. Journal Infection Disease, 2010. Vol: 202 (pp: 1212-1218).
- Waqas Nasir. *A study of Norovirus-HBGA interactions*. Chalmers University of Technology. University of Gothenburg, Sweden. December 2009.
- Arias Castellanos I and González Quiceno M. *Determinación de la coinfección entre enterobacterias, parásitos y Norovirus en niños de edades entre tres meses y cinco años con enfermedad diarreica aguda en los barrios arabia y Jerusalén de la localidad de ciudad Bolívar*. Universidad Pontificia Javeriana. Facultad de Ciencias, Bogotá (<http://www.javeriana.edu.co/biblos/tesis/ciencias/tesis372.pdf>).
- Hoonmo L. Koo and Herbert L. DuPont. *Noroviruses as a potential cause of protracted and lethal disease in immunocompromised patients*. Clinical Infectious Disease, 2009. Vol: 49 (pp: 1069-71).

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