***S. pneumoniae-L. pneumophila* MonlabTest®**

MO-804037 20 TESTS

One step test to detect *S. pneumoniae* and *L. pneumophila*

A rapid, one step test for the qualitative detection of *Streptococcus pneumoniae* and *Legionella pneumophila* serogroup 1 antigens in urine specimens.
For *in vitro* diagnostic use only.

INTENDED USE

An *in vitro* rapid immunochromatographic assay for the qualitative detection of *Streptococcus pneumoniae* and *Legionella pneumophila* serogroup 1 antigen (*L. pneumophila* serogroup 1) in urine specimens from patients with symptoms of pneumonia. It is intended to aid in the presumptive diagnosis of pneumoniae caused by *Streptococcus pneumoniae* and of *Legionella* infection (Legionnaires' Disease) caused by *L. pneumophila* serogroup 1 in conjunction with culture and other methods.

SYNTHESIS

Infection with *Streptococcus pneumoniae* is one of the leading bacterial causes of morbidity and mortality worldwide and is the principal causative agent of community-acquired pneumonia, the sixth most common cause of overall mortality in Western countries. The incidence of pneumococcal disease is highest in infants under 2 years of age and in people over 60 years of age. Recent estimates of child deaths caused by *S. pneumoniae* range from 700000 to 1 million every year worldwide. On the basis of differences in capsular polysaccharide structure *S. pneumoniae* can be divided into more than 80 serotypes that can lead to varying clinical and epidemiologic profiles.

Legionnaires' disease is a serious form of pneumonia that carries with it a mortality rate in the order of 10-15% in otherwise healthy individuals. Symptoms include a flu-like illness, followed by a dry cough and frequently progress to pneumonia. Approximately 33% of people infected may also present with diarrhoea and vomiting and around 50% may show signs of mental confusion. The incubation period normally ranges from 2-10 days with 3-6 days the typical illness onset time after exposure. Legionnaires' disease may present as an outbreak of two or more cases following a limited temporal and spatial exposure to a single source, as a series of independent cases in an area in which it is highly endemic or as sporadic cases without any obvious temporal or geographical grouping. Outbreaks have occurred repeatedly in buildings such as hotels and hospitals.

The *S. pneumoniae-L. pneumophila* MonlabTest® allows for early diagnosis of *Legionella pneumophila* serogroup 1 infection through detection of a specific soluble antigen present in the urine of patients with Legionnaires' Disease. *Legionella pneumophila* serogroup 1 antigen has been detected in urine as early as three days after the onset of symptoms. The test is rapid, giving a result within 15 minutes, and utilizes a urine specimen which is convenient for collection, transport, and subsequent detection of early, as well as later, stages of disease.

PRINCIPLE

The *S. pneumoniae-L. pneumophila* MonlabTest® is a qualitative lateral flow immunoassay for the detection of *Streptococcus pneumoniae* (*S. pneumoniae*) and *Legionella pneumophila* serogroup 1 (*L. pneumophila*) in urine samples. The membrane of the Test A is pre-coated with monoclonal antibodies against *S. pneumoniae* antigens and the membrane of the Test B is pre-coated with polyclonal antibodies against *L. pneumophila* antigens. During testing, the sample reacts with the particle coated with anti-*S. pneumoniae* antibodies in the Test A and/or with anti-*L. pneumophila* antibodies in the Test B which were pre-dried on the test strips. The mixture moves upward on the membrane by capillary action. In the case of a positive result in the Test A the specific antibodies present on the membrane will react with the mixture conjugate and generate one red colored line. In the case of a positive result in the Test B the specific antibodies present on the membrane will react with the conjugate and generate one red colored line. A blue colored line always appears in the control lines and serves as verification that sufficient volume was added, that proper flow was obtained and as an internal control for the reagents.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not use components past its expiration date.
- Do not mix components from different kit lots.
- Solutions used to make the control swabs are inactivated using standard methods. However, patient samples, controls, and Device tests should be handled as though they could transmit disease. Observe established precautions against microbial hazards.
- The test should be discarded in a proper biohazard container after testing.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

MATERIALS PROVIDED

- 20 *S. pneumoniae-L. pneumophila* MonlabTest® Tests
- 20 Testing tubes
- 20 Pipettes
- 1 Reagent (sample and controls diluent)
- Instructions for use
- 1 Positive Control Swab: Inactivated *S. pneumoniae* swab + testing tube+ pipette
- 1 Positive Control Swab: Inactivated *L. pneumophila* swab + testing tube+ pipette
- 1 Negative Control Swab: negative swab + testing tube + pipette

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

SPECIMEN COLLECTION AND PREPARATION

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30°C/59-86°F) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8°C for up to 14 days or at -10°C to -20°C for longer periods before testing.

When necessary, urine specimens should be shipped in leak-proof containers at 2-8°C or frozen.

Allow all specimens to equilibrate to room temperature before testing.


ASSAY PROCEDURE

Allow the tests, urine samples and buffer to reach room temperature (15-30°C/59-86°F) prior to testing. Do not open pouches until ready to perform the assay.

Procedure for Patient Samples. See illustration 1:

1. Use a separate testing tube for each sample. Add 6 drops of urine sample (1) and 2 drops of Reagent (2) into the testing tube and mix with the pipette.
2. Remove the device test from its pack just before use. Use a separate test for each sample. Place the test on a flat surface.
3. Dispense exactly 3 drops of the sample (urine) + Reagent, into the circular window marked with the letter A (3) and 3 drops, using the same testing tube, in the circular window marked with the letter B (4). Read the result at 15 minutes.


Illustration 1

Procedure for Positive Control Swabs. See illustration 2:

- 1 Add 13 free falling drops of Reagent into the testing tube. Hold Reagent bottle vertically (1a).
 - 2 Immediately remove *S. pneumoniae* Positive Control swab from the pouch and put the swab into the testing tube with reagent, mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab (2a).
 - 3 Remove the device from its sealed bag just before using.
 - 4 Place the test on a flat surface. Dispense exactly 3 drops from the testing tube, into the circular window marked with the letter A (3a). Read the result at 15 minutes.
- Repeat the procedure for Positive Control Swabs using *Legionella* Positive Control swab (1b, 2b), to dispense 3 drops from the testing tube, into the circular window marked with the letter B (3b).


Illustration 2

Procedure for Negative Control Swab. See illustration 3:

- 1 Add 13 free falling drops of Reagent into the testing tube. Hold Reagent bottle vertically (1).
- 2 Immediately remove Negative Control swab from the pouch and put the swab into the testing tube with reagent, mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab (2).
- 3 Remove the device from its sealed bag just before using.
- 4 Place the test on a flat surface. Dispense exactly 3 drops from the testing tube, into the circular window marked with the letter A (3) and other 3 drops from the same testing tube into the circular window marked with the letter B (4). Read the result at 15 minutes.

Illustration 3

INTERPRETATION OF RESULTS

Illustration 4

	<i>S. pneumoniae</i> and <i>L. pneumophila</i> negative.		<i>S. pneumoniae</i> positive. <i>L. pneumophila</i> negative.
	<i>S. pneumoniae</i> and <i>L. pneumophila</i> positive.		<i>S. pneumoniae</i> negative. <i>L. pneumophila</i> positive.

INVALID: Total absence of the blue control line in one or both Tests (A/B) regardless the appearance or not of the red test lines in one or both Tests (A/B). Notes: insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay using a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor. See illustration 4.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result region (T) will vary depending on the concentration in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal Quality Control:

The *S. pneumoniae-L. pneumophila* MonlabTest® contains built-in control features. The blue line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

Positive Procedural Control: The blue line at the Control line region can be considered an internal positive procedural control. If capillary flow has occurred, this line will always appear.

Negative Procedural Control: The clearing of background colour in the result window provides a negative background control.

External Positive Controls: Good Laboratory Practice recommends the use of positive controls to assure functionality of reagents and proper performance of assay procedure. Positive control swabs that will monitor the entire assay are provided in the kit.

LIMITATIONS

1. *S. pneumoniae-L. pneumophila* MonlabTest® has been validated using urine samples only. Other samples (e.g. plasma, serum or other body fluids) that may contain *Legionella* and *Streptococcus pneumoniae* antigens have not been evaluated. The test cannot be used on environmental samples (i.e. potable water).
2. A negative antigen result does not exclude infection with *S. pneumoniae* and/or *L. pneumophila* serogroup 1. Culture or other methods are recommended for suspected pneumonia to detect causative agents other than *S. pneumoniae* and *L. pneumophila* serogroup 1 and to recover *S. pneumoniae* and *L. pneumophila* serogroup 1 when antigen is not detected in urine.
3. The diagnosis of Legionnaires' disease and *Streptococcus pneumoniae* infection cannot be based on clinical or radiological evidence alone. There is no single satisfactory laboratory test for Legionnaire's and pneumococcal disease. Therefore, culture results, serology and antigen detection methods should be used in conjunction with clinical findings to make an accurate diagnosis.
4. Excretion of *Legionella* and *Streptococcus pneumoniae* antigen in urine may vary depending on the individual patient. Antigen excretion may begin as early as 3 days after onset of symptoms and persist for up to 1 year afterwards. A positive *S. pneumoniae-L. pneumophila* MonlabTest® result can occur due to current or past infection and therefore is not definitive for infection without other supporting evidence.
5. Although the test is compatible with boric acid, its presence increases the reactivity and can lead to false positives, so it is not recommended as a preservative.

EXPECTED VALUES

Worldwide each year, there are over 14 million serious *Streptococcus pneumoniae* infections in children <5 years of age leading to over 800,000 deaths. The disease occurs in both epidemic and endemic forms and sporadic cases are not easily differentiated from other respiratory infections by clinical symptoms. An estimated 25,000 to 100,000 cases of *Legionella* infection occur in the United States annually. The resulting mortality rate, ranging from 25% to 40%, can be lowered if the disease is diagnosed rapidly and appropriate antimicrobial therapy is instituted early.

PERFORMANCE CHARACTERISTICS

Sensitivity and specificity (Test A)

The *S. pneumoniae-L. pneumophila* MonlabTest® was highly specific (>98%) and also sensitive (>89%) compared with the results of the Binax NOW® *Streptococcus pneumoniae* Antigen Card.

Sensitivity and specificity (Test B)

The *S. pneumoniae-L. pneumophila* MonlabTest® was highly specific (>99%) and also sensitive (>99%) compared with the results of the Binax NOW® *Legionella* Urinary Antigen Test.

The *S. pneumoniae-L. pneumophila* MonlabTest® was used to evaluate 111 samples for strip A evaluation (*S. pneumoniae*) vs BinaxNOW® *Streptococcus pneumoniae* Antigen Card, Binax.

The *S. pneumoniae-L. pneumophila* MonlabTest® was used to evaluate 149 samples for strip B evaluation (*Legionella*) vs. BinaxNOW® *Legionella* Urinary Antigen Card, Binax.

Cross reactivity

It was performed an evaluation to determine the cross reactivity of *S. pneumoniae-L. pneumophila* MonlabTest®. There is no cross reactivity with other pathogens occasionally present in urine:

- *Legionella pneumophila* (Strip A)
- *S. pneumoniae* (Strip B)

REFERENCES

1. Roig, J., X. Aquiler, J. Ruiz, et al. Comparative study of *Legionella pneumophila* and other nosocomial-acquired pneumonias. Chest. 1991;99:344-50.
2. Berdal, B.P., C.E. Farshy, and J.C. Feeley. Detection of *Legionella pneumophila* antigen in urine by enzyme-linked immunospecific assay. J. Clin. Microbiol. 1979;9:575-578.
3. White A., et al. Rapad diagnosis of Legionnaires' disease. Trans Am Clin. Climatol. Assoc. 1982;93:50-62
4. Bibb, W.F., P.M. Arnow, L. Thacker, and R.M. McKinney. Detection of soluble *Legionella pneumophila* antigens in serum and urine specimens by enzyme-linked immunosorbent assay with monoclonal and polyclonal antibodies. J. Clin. Microbiol. 1984;20:478-482.
5. Tang, P.W., and S. Toma. Broad-spectrum enzyme-linked immunosorbent assay for detection of *Legionella* soluble antigens. J. Clin. Microbiol. 1986;24:556-558.
6. Kohler, R.B., W.C. Winn, Jr., and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984;20:605-607.
7. Nationwide Trends in Bacterial Meningitis before the Introduction of 13-Valent Pneumococcal Conjugate Vaccine-Burkina Faso, 2011-2013. PLoS One. 2016 Nov 10; 11(11).
8. Burden of disease caused by *Streptococcus pneumoniae* in children younger than 5 years: global estimates. Lancet. 2009 Sep 12;374(9693):893-902.
9. MIERNYK, KAREN M.; BULKOW, LISA R.; CASE, SAMANTHA L.; ET AL. "Population structure of invasive *Streptococcus pneumoniae* isolates among Alaskan children in the conjugate vaccine era, 2001 to 2013". DIAGNOSTIC MICROBIOLOGY AND INFECTIOUS DISEASE 2016; 86(2): 224-230.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only		Positive control swab
	Don't re-use		Consult instructions for use		Negative control swab
	Contains sufficient for <n> tests		Keep dry		
	Catalogue Code		Temperature limitation		
	Lot Number		Use by		

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- *Legionella pneumophila* (Strip A)

- *S. pneumoniae* (Strip B)

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1. Roig, J., X. Aquiler, J. Ruiz, et al. Comparative study of *Legionella pneumophila* and other nosocomial-acquired pneumonias. Chest. 1991;99:344-50.
2. Berdal, B.P., C.E. Farshy, and J.C. Feeley. Detection of *Legionella pneumophila* antigen in urine by enzyme-linked immunospecific assay. J. Clin. Microbiol. 1979;9:575-578.
3. White A., et al. Rapid diagnosis of Legionnaires' disease. Trans Am Clin. Climatol. Assoc. 1982;93:50-62
4. Bibb, W.F., P.M. Arnow, L. Thacker, and R.M. McKinney. Detection of soluble *Legionella pneumophila* antigens in serum and urine specimens by enzyme-linked immunosorbent assay with monoclonal and polyclonal antibodies. J. Clin. Microbiol. 1984;20:478-482.
5. Tang, P.W., and S. Toma. Broad-spectrum enzyme-linked immunosorbent assay for detection of *Legionella* soluble antigens. J. Clin. Microbiol. 1986;24:556-558.
6. Kohler, R.B., W.C. Winn, Jr., and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984;20:605-607.
7. Nationwide Trends in Bacterial Meningitis before the Introduction of 13-Valent Pneumococcal Conjugate Vaccine-Burkina Faso, 2011-2013. PLoS One. 2016 Nov 10; 11(11).
8. Burden of disease caused by *Streptococcus pneumoniae* in children younger than 5 years: global estimates. Lancet. 2009 Sep 12;374(9693):893-902.
9. MIERNYK, KAREN M.; BULKOW, LISA R.; CASE, SAMANTHA L.; ET AL. "Population structure of invasive *Streptococcus pneumoniae* isolates among Alaskan children in the conjugate vaccine era, 2001 to 2013". DIAGNOSTIC MICROBIOLOGY AND INFECTIOUS DISEASE 2016; 86(2): 224-230.

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