**S. Pneumoniae MonlabTest®**

**MO-804036 20 TESTS**

*One step test to detect Streptococcus pneumoniae*

A rapid, one step test for the qualitative detection of *Streptococcus pneumoniae* in urine specimens.

For *in vitro* diagnostic use only.

**INTENDED USE**

An *in vitro* rapid immunochromatographic assay for the qualitative detection of *Streptococcus pneumoniae* in urine specimens from patients with symptoms of pneumonia. It is intended to aid in the presumptive diagnosis of pneumoniae caused by *Streptococcus pneumoniae* in conjunction with culture and other methods.

**SYNTHESIS**

Infection with *Streptococcus pneumoniae* is one of the leading bacterial causes of morbidity and mortality worldwide and is the principal causative agent of community-acquired pneumonia, the sixth most common cause of overall mortality in Western countries. The incidence of pneumococcal disease is highest in infants under 2 years of age and in people over 60 years of age. Recent estimates of child deaths caused by *S pneumoniae* range from 700000 to 1 million every year worldwide. On the basis of differences in capsular polysaccharide structure *S pneumoniae* can be divided into more than 80 serotypes that can lead to varying clinical and epidemiologic profiles.

**PRINCIPLE**

The *S. Pneumoniae* MonlabTest® is an immunochromatographic membrane assay to detect *Streptococcus pneumoniae* antigen in human urine sample. Anti-*Streptococcus pneumoniae* antibody, the test line, is adsorbed onto nitrocellulose membrane. Antibodies of the control line were adsorbed onto the same membrane as a second band. Anti-*Streptococcus pneumoniae* antibodies are conjugated to visualizing particles that are dried onto an inert absorbent support.

During testing the sample is allowed to react with the conjugate which was pre-adsorbed on the strip test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. *S. pneumoniae* antigen captured by immobilized anti-*S. pneumoniae* antibody reacts to bind conjugated antibody. The other immobilized antibodies also captures visualizing conjugate, forming the control line. A positive test result is read visually in 15 minutes or less depending on the concentration of antigen present in the urine specimen. A negative *S. pneumoniae* result, read in 15 minutes, indicates that *S. pneumoniae* antigen was not detected in the urine sample.

The test is interpreted by the presence or absence of visually detectable red and blue colored lines. A positive result will include the detection of both a test and a control line, while a negative assay will produce only the control line. Failure of the control line to appear, whether the test line is present or not, indicates an invalid assay.

**PRECAUTIONS**

- For professional *in vitro* diagnostic use only.
- Do not use components past its expiration date.
- Do not mix components from different kit lots.
- Solutions used to make the control swabs are inactivated using standard methods. However, patient samples, controls, and Device tests should be handled as though they could transmit disease. Observe established precautions against microbial hazards.
- The test should be discarded in a proper biohazard container after testing.

**STORAGE AND STABILITY**

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
<ul style="list-style-type: none"> <li>- 20 <i>S. Pneumoniae</i> MonlabTest® Tests</li> <li>- Reagent (sample and controls diluent)</li> <li>- Instructions for use</li> <li>- 1 Positive Control Swab: Inactivated <i>S. pneumoniae</i> swab + testing tube+ pipette</li> <li>- 1 Negative Control Swab: <i>S. pneumoniae</i> negative swab + testing tube + pipette</li> <li>- 20 Testing tubes</li> <li>- 20 Plastic pipettes</li> </ul>	<ul style="list-style-type: none"> <li>- Specimen collection container</li> <li>- Disposable gloves</li> <li>- Timer</li> </ul>

**SPECIMEN COLLECTION AND PREPARATION**

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30°C/59-86°F) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8°C for up to 14 days or at -10°C to -20°C for longer periods before testing.

When necessary, urine specimens should be shipped in leak-proof containers at 2-8°C or frozen.

Allow all specimens to equilibrate to room temperature before testing.

**ASSAY PROCEDURE**

**Procedure for Patient Samples. See illustration 1:**


Do not remove device test from pack until test sample has reached room temperature. Bring patient urine and/or liquid urine control(s) to room temperature (15-30°C/59-86°F).

1. Use a separate testing tube for each sample. Add 3 drops of urine sample and 1 drop of Reagent into the testing tube and mix with the pipette.
2. Remove the device test from its pack just before use. Use a separate test for each sample. Place the test on a flat surface
3. Dispense exactly 3 drops from the testing tube, into the circular window marked with S. Start the timer. Read the result at **15 minutes**.


**Procedure for Positive and Negative Swab Controls. See illustration 2:**

1. Add slowly 13 drops of Reagent into the testing tube.
  2. Immediately remove swab control from the pouch and put the swab into the testing tube with reagent, mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab.
  3. Remove the device test from its sealed bag just before using and place the test on a flat surface. Dispense exactly 3 drops from the testing tube, into the circular window marked with S. Start the timer. Read the result at **15 minutes**.
- Repeat the procedure for Negative Swab Control using the same Reagent used for test and for positive control swab.

**Illustration 1**


**Illustration 2**

### INTERPRETATION OF RESULTS

#### Illustration 3

**POSITIVE** sample or positive control swab: two lines appear across the central window, in the result line region, **red** test line marked with the letter T, and in the control line region, **blue** control line marked with the letter C. This means that antigen was detected.

**NEGATIVE** sample or negative control swab: will give only the **blue** band across the control line region marked with the letter C (control line). The control line means that the detection part of the test was done correctly, but no *S. pneumoniae* antigen was detected.

**INVALID**: if no lines are seen, or if just the test line is seen, the assay is invalid. Invalid tests should be repeated. If the problem persists, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result region (T) will vary depending on the concentration in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

### QUALITY CONTROL

#### Internal Quality Control:

The *S. pneumoniae* MonlabTest® contains built-in control features. The blue line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

**Positive Procedural Control:** The blue line at the Control line region can be considered an internal positive procedural control. If capillary flow has occurred, this line will always appear.

**Negative Procedural Control:** The clearing of background colour in the result window provides a negative background control.

**External Positive and Negative Controls:** Good Laboratory Practice recommends the use of positive and negative controls to assure functionality of reagents and proper performance of assay procedure. Positive and negative control swabs that will monitor the entire assay are provided in the kit.

### LIMITATIONS

- S. pneumoniae* MonlabTest® has been validated using human urine samples only. Other samples (e.g. plasma, serum or other body fluids) that may contain *Streptococcus pneumoniae* antigen have not been evaluated. The test cannot be used on environmental samples (i.e. potable water).
- A negative antigen result does not exclude infection with *S. pneumoniae*. Culture or other methods are recommended for suspected pneumoniae to detect causative agents other than *S. pneumoniae* and to recover *S. pneumoniae* when antigen is not detected in urine.
- The diagnosis of *Streptococcus pneumoniae*'s disease cannot be based on clinical or radiological evidence alone. There is no single satisfactory laboratory test for *Streptococcus pneumoniae*'s disease. Therefore, culture results, serology and antigen detection methods should be used in conjunction with clinical findings to make an accurate diagnosis.
- Excretion of *Streptococcus pneumoniae* antigen in urine may vary depending on the individual patient. A positive *S. pneumoniae* MonlabTest® result can occur due to current or past infection and therefore is not definitive for infection without other supporting evidence.

### EXPECTED VALUES

Worldwide each year, there are over 14 million serious Streptococcus pneumoniae infections in children <5 years of age leading to over 800,000 deaths.

### PERFORMANCE CHARACTERISTICS

#### Sensitivity and specificity

The *S. pneumoniae* MonlabTest® was highly specific (>98%) and also sensitive (>89%) compared with the results of the Binax NOW® *Streptococcus pneumoniae* Antigen Card.

Sensitivity >89% and Specificity >98%

#### Cross reactivity



It was performed an evaluation to determine the cross reactivity of *S. pneumoniae* MonlabTest®. There is no cross reactivity with other pathogens occasionally present in urine:

- *Legionella pneumophila*

### REFERENCES

- Nationwide Trends in Bacterial Meningitis before the Introduction of 13-Valent Pneumococcal Conjugate Vaccine-Burkina Faso, 2011-2013. PLoS One. 2016 Nov 10; 11(11).
- Burden of disease caused by Streptococcus pneumoniae in children younger than 5 years: global estimates. Lancet. 2009 Sep 12;374(9693):893-902.
- MIERNYK, KAREN M.; BULKOW, LISA R.; CASE, SAMANTHA L.; ET AL. "Population structure of invasive Streptococcus pneumoniae isolates among Alaskan children in the conjugate vaccine era, 2001 to 2013". DIAGNOSTIC MICROBIOLOGY AND INFECTIOUS DISEASE 2016; 86(2): 224-230.

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only		Positive control swab
	Don't re-use		Consult instructions for use		Negative control swab
	Contains sufficient for <n> tests		Keep dry		
	Catalogue Code		Temperature limitation		
	Lot Number		Use by		

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- *Legionella pneumophila*

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- Nationwide Trends in Bacterial Meningitis before the Introduction of 13-Valent Pneumococcal Conjugate Vaccine-Burkina Faso, 2011-2013. PLoS One. 2016 Nov 10; 11(11).
- Burden of disease caused by *Streptococcus pneumoniae* in children younger than 5 years: global estimates. Lancet. 2009 Sep 12;374(9693):893-902.
- MIERNYK, KAREN M.; BULKOW, LISA R.; CASE, SAMANTHA L.; ET AL. "Population structure of invasive *Streptococcus pneumoniae* isolates among Alaskan children in the conjugate vaccine era, 2001 to 2013". DIAGNOSTIC MICROBIOLOGY AND INFECTIOUS DISEASE 2016; 86(2): 224-230.

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